



MedRAX Conflict Resolution: Complete Pipeline Analysis

High-Level Architecture








Complete Example Scenario: Patient with Suspected Cardiomegaly

Stage 1: Tool Outputs (Raw)



REPORT GENERATOR (ViT-BERT)

Output: "CHEST X-RAY REPORT

FINDINGS:

The cardiac silhouette is mildly enlarged. Lungs are clear.

IMPRESSION:

Mild cardiomegaly. No acute cardiopulmonary process."

Metadata: {overall_confidence: 0.82}

CHEXAGENT VQA (CheXagent-2-3B)

Question: "Is there cardiomegaly in this image?"

Output: "The heart size appears within normal limits. No evidence of cardiomegaly is seen on this frontal chest radiograph."

Metadata: {self_consistency_score: 0.68}

CLASSIFIER (DenseNet-121)

```
Output: {
  "Atelectasis": 0.12,
  "Cardiomegaly": 0.22,      ← 22% probability
  "Consolidation": 0.08,
  "Edema": 0.05,
  "Effusion": 0.15,
  "Pneumonia": 0.11,
  "Pneumothorax": 0.03,
  ... (18 pathologies total)
}
```

SEGMENTATION (PSPNet)

```
Output: {
  "segmentation_image_path": "temp/segmentation_8022ced0.png",
  "metrics": {
    "Heart": {
      "area_pixels": 58000,
```

```
| | "area_cm2": 23.2,
| |
| | "width": 185,           ← Heart width in pixels
| |
| | "height": 168,
| |
| | "confidence_score": 0.94
| |
| | },
| |
| | "Left Lung": {"width": 162, ...},
| |
| | "Right Lung": {"width": 178, ...}
| |
| | }
| |
| | }
```

Stage 2: Normalization to CanonicalFinding

STAGE 2: NORMALIZED FINDINGS

All tool outputs converted to uniform CanonicalFinding format:

```
| Finding 1: Report Generator
|
| pathology: "Cardiomegaly"
|
| source_tool: "chest_xray_report_generator"
|
| confidence: 0.82
```

```
| | evidence_type: "report"
| | raw_value: "Mild cardiomegaly..."
| | metadata: {text: "cardiac silhouette is mildly enlarged..."}
```

```
| Finding 2: CheXagent VQA
| | pathology: "Cardiomegaly"
| | source_tool: "chest_xray_expert"
| | confidence: 0.25    ← Low because says "normal"
| | evidence_type: "vqa"
| | raw_value: "No evidence of cardiomegaly..."
| | metadata: {text: "heart size appears within normal limits..."}
```

```
| Finding 3: Classifier
| | pathology: "Cardiomegaly"
| | source_tool: "chest_xray_classifier"
| | confidence: 0.22
| | evidence_type: "classification"
| | raw_value: {"Cardiomegaly": 0.22, ...}
| | metadata: {text: "Cardiomegaly probability: 22%"}
```

```
| Finding 4: Segmentation
|
|   pathology: "Cardiomegaly"
|
|   source_tool: "chest_xray_segmentation"
|
|   confidence: 0.89    ← High because CTR > 0.5
|
|   evidence_type: "segmentation"
|
|   raw_value: {"Heart": {"width": 185}, "Left Lung": {"width": 162}...}
|
|   metadata: {cardiothoracic_ratio: 0.544}
```

Stage 3: Conflict Detection (Three-Step Process)

STAGE 3: CONFLICT DETECTION

STEP 1: BERT (Text vs Text)

Comparing: Report Generator ↔ CheXagent VQA

Text 1: "cardiac silhouette is mildly enlarged...Mild cardiomegaly"

Text 2: "heart size appears within normal limits. No evidence of cardiomegaly"|

BERT NLI Output:

	contradiction_prob: 0.91	✓ HIGH
	entailment_prob: 0.03	
	neutral_prob: 0.06	

Result: ✓ CONFLICT #1 DETECTED

Type: semantic

Severity: critical (0.91 > 0.85)

Tools: [report_generator, chest_xray_expert]

STEP 2: Rule-Based (Confidence Gap)

Confidences: [0.82, 0.25, 0.22, 0.89]

Max: 0.89 (Segmentation)

Min: 0.22 (Classifier)

Gap: 0.89 - 0.22 = 0.67 > 0.4 threshold ✓

Check: Max (0.89) > 0.7? YES

Check: Min (0.22) < 0.3? YES

Result: ✓ CONFLICT #2 DETECTED (but deduplicated - BERT caught similar)

STEP 3: GACL (Measurements vs Probabilities)

Segmentation Measurements:

Heart width: 185 pixels
Thorax width: 162 + 178 = 340 pixels
CTR: 185 / 340 = 0.544

Medical Threshold: CTR > 0.50 = Cardiomegaly

GACL Analysis:

	Measurements say: CARDIOMEGALY (CTR = 54.4% > 50%)	
	Classifier says: NO CARDIOMEGALY (22% probability)	
	INCONSISTENCY DETECTED!	

Result: ✓ CONFLICT #3 DETECTED

Type: semantic (anatomical)
Tools: [segmentation_tool, classifier]

Explanation: "CTR of 0.544 suggests cardiomegaly but classifier reports 22%"

Stage 4: Conflict Resolution

STAGE 4: CONFLICT RESOLUTION

RESOLVING CONFLICT #1: Report Generator vs CheXagent

Step 4.1: Analyze BERT Scores

contradiction_prob: 0.91 (HIGH)

entailment_prob: 0.03 (LOW)

is_false_positive: NO (entailment < 0.7)

severity_adjustment: 1.0 (no discount needed)

Step 4.2: BERT-Guided Resolution

- BERT contradiction > 0.85? YES (0.91)
- Confidence gap: 0.82 - 0.25 = 0.57 > 0.3? YES

- Clear winner: Report Generator (0.82)

Resolution #1:

	decision: "bert_high_confidence_leader"	
	selected_tool: "chest_xray_report_generator"	
	value: TRUE (Cardiomegaly present)	
	confidence: 0.82	
	should_defer: FALSE	

RESOLVING CONFLICT #3: Segmentation vs Classifier (GACL)

Step 4.3: Task-Aware Arbitration

- Conflict type: semantic (anatomical)
- For anatomical measurements: Trust segmentation
- Segmentation confidence: 0.89
- Classifier confidence: 0.22

Resolution #3:

	decision: "trust_segmentation_measurements"	
	selected_tool: "chest_xray_segmentation"	
	value: TRUE (Cardiomegaly present)	

```
| confidence: 0.89 |
| reasoning: "CTR=0.544 objectively exceeds 0.50" |
| should_defer: FALSE |
└──────────────────┘
```

Stage 5: Final Output

STAGE 5: FINAL OUTPUT

⚠ CONFLICT DETECTION REPORT

Detected 2 conflict(s)

Timestamp: 2026-02-05 14:30:00

Conflict #1 - CRITICAL SEVERITY

Type: semantic

Finding: Cardiomegaly

Tools: chest_xray_report_generator, chest_xray_expert

• report_generator: "Mild cardiomegaly" (confidence: 82.0%)

• chest_xray_expert: "No evidence of cardiomegaly" (confidence: 25.0%)

Resolution:

Decision: bert_high_confidence_leader

Selected: chest_xray_report_generator

Confidence: 82.0%

Reasoning: BERT detected high contradiction (91%). Trusting report generator with significantly higher confidence (82% vs 25%).

Conflict #2 - CRITICAL SEVERITY

Type: semantic (anatomical)

Finding: Anatomical pattern consistency

Tools: segmentation_tool, chest_xray_classifier

- segmentation: CTR=0.544 (confidence: 89.0%)
- classifier: Cardiomegaly=22% (confidence: 22.0%)

Resolution:

Decision: trust_segmentation_measurements

Selected: chest_xray_segmentation

Confidence: 89.0%

Reasoning: Graph-based anatomical analysis shows CTR=0.544 > 0.50 threshold. Objective measurements override classifier.

FINAL DIAGNOSIS: CARDIOMEGALY PRESENT

Confidence: 85.5% (average of 82% + 89%)

Supporting Evidence:

✓ Report Generator: "Mild cardiomegaly" (82%)

✓ Segmentation: CTR = 0.544 > 0.50 (89%)

Contradicting Evidence:

✗ CheXagent VQA: "No cardiomegaly" (25%)

✗ Classifier: 22% probability

Status: ✓ RESOLVED (No human review needed)

Summary: What Each Component Does

COMPONENT RESPONSIBILITIES	
BERT (DeBERTa-MNLI)	<div><div>└─ Input: Two text strings</div><div>└─ Output: contradiction/entailment/neutral probabilities</div><div>└─ Detects: Text-based contradictions</div><div>└─ Example: "enlarged heart" vs "normal heart" → 91% contradiction</div></div>
GACL (Graph-Based Anatomical Consistency)	<div><div>└─ Input: Segmentation metrics + Classifier probabilities</div><div>└─ Output: Conflict detected or not + explanation</div></div>

- └─ Detects: Measurement vs probability inconsistencies

└─ Example: CTR=0.544 vs Cardiomegaly=22% → CONFLICT
- Rule-Based
- └─ Input: List of confidence scores

└─ Output: Presence conflict detected or not

└─ Detects: Large confidence gaps (one says YES, another says NO)

└─ Example: 89% vs 22% gap = 67% > 40% threshold → CONFLICT
- ConflictResolver
- └─ Input: Detected conflicts + all findings


└─ Output: Resolution decision with reasoning

└─ Methods: BERT scores → Tool expertise → Weighted average → Defer

└─ Example: High BERT contradiction + confidence leader → Trust winner

Coverage Matrix

WHAT EACH DETECTOR CAN COMPARE		
	BERT	GACL
Report ↔ CheXagent	✓	✗
Report ↔ LLaVA-Med	✓	✗
CheXagent ↔ LLaVA	✓	✗



Models Used

Tool	Model	Output Type
Report Generator	<code>IAMJB/chexpert-mimic-cxr-*-baseline</code> (ViT-BERT)	Text
CheXagent VQA	<code>StanfordAIMI/CheXagent-2-3b</code>	Text
Classifier	<code>torchxrayvision.models.DenseNet</code> (densenet121-res224-all)	18 Probabilities
Segmentation	<code>torchxrayvision.baseline_models.chestx_det.PSPNet</code>	14 Organ Metrics
LLaVA-Med	<code>microsoft/llava-med-v1.5-mistral-7b</code>	Text
BERT Conflict Detector	<code>microsoft/deberta-v3-large-mnli</code>	NLI Scores

Key Files

File	Purpose
<code>medrax/agent/conflict_resolution.py</code>	Conflict detection & resolution logic
<code>medrax/agent/bert_conflict_detector.py</code>	BERT NLI-based text comparison
<code>medrax/agent/anatomical_consistency_graph.py</code>	GACL measurement analysis
<code>medrax/agent/canonical_output.py</code>	Output normalization
<code>medrax/tools/classification.py</code>	DenseNet classifier
<code>medrax/tools/segmentation.py</code>	PSPNet segmentation
<code>medrax/tools/report_generation.py</code>	ViT-BERT report generator
<code>medrax/tools/xray_vqa.py</code>	CheXagent VQA